

Product Information

MemDX™ Membrane Protein Human ADTRP (Androgen dependent TFPI regulating protein)

Full Length

Cat. No.: **MPC4414K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ADTRP membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ADTRP

Protein Length

Full length

Protein Class

Receptor

TMD

6

Sequence

MTKTSTCIYHFLVLSWYTFLNYYISQEGKDEVKPKILANGARWKYMTLLN
LLLQTIYFGVTCDDVLKRTKGGKDIKFLTAFRDLLFTTLAFPVSTFVFL
AFWILFLYNRDLIYPKVLDTVIPVWLNHAMHTFIFPITLAEVVLRRPHSYP
SKKTGLTLLAAASIAYISRILWLYFETGTWVYPVFAKLSLLGLAAFFSLS
YVFIASIYLLGEKLNHWKWGDMRQPRKKRK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

ADTRP

Full Name

Androgen dependent TFPI regulating protein

Alternative Names

ADTRP; AIG1L; C6orf105; dJ413H6.1; androgen-dependent TFPI-regulating protein; FAHFA hydrolase ADTRP; androgen-dependent TPF1-regulating protein; fatty acid esters of hydroxy fatty acids hydrolase ADTRP; Androgen dependent TFPI regulating protein

Gene ID

[84830](#)

UniProt ID

[Q96IZ2](#)