

Product Information

MemDX™ Membrane Protein Human SMIM3 (Small integral membrane protein 3) Full Length

Cat. No.: **MPC4418K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SMIM3 membrane protein expressed in *E.coli*. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SMIM3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MDAVSQVPMEVVLPHKILDIWVIVLILATIVIMTSLLLCPATAVIYRM
RTHPILSGAV

Product Description

Expression Systems

E.coli

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw

cycles.

Target

Target Protein

SMIM3

Full Name

Small integral membrane protein 3

Alternative Names

SMIM3; NID67; MST150; C5orf62; NGF-induced differentiation clone 67 protein; Uncharacterized protein C5orf62; small membrane protein NID67; Small integral membrane protein 3

Gene ID

[85027](#)

UniProt ID

[Q9BZL3](#)