

Product Information

MemDX™ Membrane Protein Human IZUMO3 (IZUMO family member 3) Full Length

Cat. No.: **MPC4422K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human IZUMO3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

IZUMO3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MGDLWLFLLLPLSAFHGVKGCLECDPKFIEDVGSLLGNLIPSEVPGRTQL
LERQIKEMIHLSEFKVSHSDKRLRLVLAVQQVVKLRWLKNEFYKLGNETWK
GVFIYQGKLLDVCQNLESKLKELLKNFSEIACSEDCIVVEGPILDCWTCL
RMTNRCFKGEYCGDEDPRKAENREIALFLILLATAVILGSAVLLFHFCIF
HRRKMKAIRRSLKEYVEKKLEELMGKIDEKEEKDFRLRK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

IZUMO3

Full Name

IZUMO family member 3

Alternative Names

IZUMO3; C9orf134; bA20A20.1; izumo sperm-egg fusion protein 3; IZUMO family member 3

Gene ID

[100129669](#)

UniProt ID

[Q5VZ72](#)