

# Product Information

## MemDX™ Membrane Protein Human ZDHHC19 (Zinc finger DHHC-type palmitoyltransferase

### 19) Full Length

Cat. No.: **MPC4569K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC19 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

ZDHHC19

#### Protein Length

Full length

#### Protein Class

Transferase

#### TMD

4

#### Sequence

MTLLTDATPLVKEPHPLPLVPRWFLPSLFAAFNVVLLVFFSGLFFAFPC  
RWLAQNGEWAFFVITGSLFVLTFVSLVSLNFSDPGILHQGSAEQGPLTVH  
VWVNHGAFRLQWCPKCCFHRPPRTYHCPWCNICVEDFDHHCKWVNNCIG  
HRNFRFFMLLVLSCLYSGAMLVTCILFLVRTTHLPFSTDKAIAIVVAVS  
AAGLLVPLSLLLLIQALSVSSADRITYKGKCRHLQGYNPFDQGCASNWYLT  
ICAPLGPKYMAEAVQLQRVVGPDWTSMPNLHPPMSPSALNPPAPTSGSLQ  
SREGTPGAW

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

ZDHHC19

**Full Name**

Zinc finger DHHC-type palmitoyltransferase 19

**Alternative Names**

ZDHHC19; DHHC19; palmitoyltransferase ZDHHC19; DHHC-19; probable palmitoyltransferase ZDHHC19; zinc finger DHHC domain-containing protein 19; zinc finger DHHC-type containing 19; zinc finger, DHHC domain containing 19; Zinc finger DHHC-type palmitoyltransferase 19

**Gene ID**

[131540](#)

**UniProt ID**

[Q8WVZ1](#)