

Product Information

MemDX™ Membrane Protein Human ZDHHC22 (Zinc finger DHHC-type palmitoyltransferase

22) Full Length

Cat. No.: **MPC4572K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC22 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC22

Protein Length

Full length

Protein Class

Transferase

TMD

5

Sequence

MLALRLN VVAPAYFLCISLVTFVLQLFLFLPSMREDPAAARLFSPALLH
GALFLFLSANALGNYVLVIQNSPDDL GACQGASARKTPCSPSTHFCRVC
ARVTLRHDHHCFFTGNCIGSRNMRNFVLFCLYTSACLYSMVAGVAYISA
VLSISFAHPLAFLTLLPTSISQFFSGAVLGSEMFVILMLYLWFAIGLACA
GFCCHQLLLILRGQTRHQVRKGVAVRARPWKRNQLQEVFGKRWLLGLLVPM
FNVGSESSKQQDK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC22

Full Name

Zinc finger DHHC-type palmitoyltransferase 22

Alternative Names

ZDHHC22; C14orf59; palmitoyltransferase ZDHHC22; DHHC-22; putative palmitoyltransferase ZDHHC22; zinc finger DHHC domain-containing protein 22; zinc finger DHHC-type containing 22; zinc finger, DHHC domain containing 22; Zinc finger DHHC-type palmitoyltransferase 22

Gene ID

[283576](#)

UniProt ID

[Q8N966](#)