

Product Information

MemDX™ Membrane Protein Human ZDHHC12 (Zinc finger DHHC-type palmitoyltransferase

12) Full Length

Cat. No.: **MPC4574K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC12 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC12

Protein Length

Full length

Protein Class

Transferase

TMD

4

Sequence

MAPWALLSPGVLVRTGHTVLTWGITLVFLHDTLRQWEEQGELLPLTF
LLLVLGSLLLYLAVSLMDPGYVNVQPQPQEELKEEQTAMVPPAIPLRRCR
YCLVLQPLRARHCRECRRCVRRYDHHCPWMENCVGERNHPLFVVYLALQL
VLLWGLYLAWSGLRFFQPWGQWLRSSGLLFATFLLLSLFSLVASLLLV
HLYLVASNTTWEFISSHRIAYLRQRPSNPFDRGLTRNLAHFFCGWPSGS
WETLWAEIEEEEGSSPAV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC12

Full Name

Zinc finger DHHC-type palmitoyltransferase 12

Alternative Names

ZDHHC12; ZNF400; DHHC-12; palmitoyltransferase ZDHHC12; DHHC domain-containing cysteine-rich protein 12; probable palmitoyltransferase ZDHHC12; zinc finger DHHC domain-containing protein 12; zinc finger DHHC-type containing 12; zinc finger protein 400; Zinc finger DHHC-type palmitoyltransferase 12

Gene ID

[84885](#)

UniProt ID

[Q96GR4](#)