

Product Information

MemDX™ Membrane Protein Human ZDHHC21 (Zinc finger DHHC-type palmitoyltransferase

21) Full Length

Cat. No.: **MPC4575K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC21 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC21

Protein Length

Full length

Protein Class

Transferase

TMD

4

Sequence

MGLRIHFVVDPHGWCCMGLIVFWLYNIVLIPKIVLFPHYEEGHIPGILI
IIFYGISIFCLVALVRASITDPGRLPENPKIPHGEREFWELCNKCNLMRP
KRSHHCSRCHGHCVRMDHHCPWINNCVGEDNHWLFLQLCFYTELLTCYAL
MFSFCHYYYYFLPLKKNLDLFVFRHELAIMRLAAFMGITMLVGITGLFYT
QLIGIITDTTSIEKMSNCCEDISRPRKPWQQTFSEVFGTRWKILWFIPFR
QRQPLRVPHYHFANHV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC21

Full Name

Zinc finger DHHC-type palmitoyltransferase 21

Alternative Names

ZDHHC21; DNZ1; DHHC21; DHHC-21; HSPC097; palmitoyltransferase ZDHHC21; DHHC domain-containing cysteine-rich protein 21; zinc finger DHHC-type containing 21; zinc finger, DHHC domain containing 21; Zinc finger DHHC-type palmitoyltransferase 21

Gene ID

[340481](#)

UniProt ID

[Q8IVQ6](#)