

Product Information

MemDX™ Membrane Protein Human ZDHHC18 (Zinc finger DHHC-type palmitoyltransferase

18) Full Length

Cat. No.: **MPC4576K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC18 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC18

Protein Length

Full length

Protein Class

Transferase

TMD

4

Sequence

MKDCEYQQISPGAAPLPASPGARRPGPAASPTPGPGPAPPAAPAPPRWSS
SGSGSGSGSGSLGRRPRRKWEVFPGRNRFYCGGRLMLAGHGGVFALTLL
ILTTTGLFFVFDCPYLARKLTLAIPHAILFFFVMSCLLQTSFTDPGIL
PRATVCEAAALEKQIDNTGSSTYRPPPRRTREVLINGQMVKLKYCFTCKMF
RPPRTSHCSVCDNCVERFDHHCPCWVGNCVGRNRYFFYAFILSLSFLTAF
IFACVVTHLTLRAQGSNFLSTLKETPASVLELVICFFSIWSILGLSGFHT
YLVASNLTTNEDIKGSWSSKRGGEASVNPYSHKSIITNCCAVLCGPLPPS
LIDRRGFVQSDTVLPSPIRSDEPACRAKPDASMVGGHP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC18

Full Name

Zinc finger DHHC-type palmitoyltransferase 18

Alternative Names

ZDHHC18; DHHC18; DHHC-18; palmitoyltransferase ZDHHC18; DHHC domain-containing cysteine-rich protein 18; zinc finger DHHC domain-containing protein 18; zinc finger DHHC-type containing 18; zinc finger, DHHC domain containing 18; Zinc finger DHHC-type palmitoyltransferase 18

Gene ID

[84243](#)

UniProt ID

[Q9NUE0](#)