

Product Information

MemDX™ Membrane Protein Human ARL8A (ADP ribosylation factor like GTPase 8A) Full Length

Cat. No.: **MPC4580K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ARL8A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ARL8A

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MIALFNKLLDWFKALFWKEEMELTLVGLQYSGKTTFVNVIASGQFNEDMI
PTVGFNMRKITKGNVTIKLWDIGGQPRFRSMWERYCRGVSAIVYMVDAAD
QEKEASKNELHNLLDKPQLQGIPVLVLGNKRDLPALDEKELIEKMNLS
AIQDREICCYSISCKEKDNIDITLQWLIQHSKSRRS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ARL8A

Full Name

ADP ribosylation factor like GTPase 8A

Alternative Names

ARL8A; GIE2; ARL10B; ADP-ribosylation factor-like protein 8A; ADP-ribosylation factor-like 10B; ADP-ribosylation factor-like 8A; novel small G protein indispensable for equal chromosome segregation 2; ADP ribosylation factor like GTPase 8A

Gene ID

[127829](#)

UniProt ID

[Q96BM9](#)