

Product Information

MemDX™ Membrane Protein Human BET1L (Bet1 golgi vesicular membrane trafficking protein like) Full Length

Cat. No.: **MPC4582K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human BET1L membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

BET1L

Protein Length

Full length

Protein Class

Transporter

TMD

1

Sequence

MADWARAQSPGAVEEILDRENKRMADSLASKVTRLKSLALDIDRDAEDQN
RYLDGMDSDFTSMTSLLTGSVKRFSTMARSGQDNRKLLCGMAVGLIVAFF
ILSYFLSRART

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

BET1L

Full Name

Bet1 golgi vesicular membrane trafficking protein like

Alternative Names

BET1L; GS15; BET1L1; GOLIM3; HSPC197; BET1-like protein; GOS-15; blocked early in transport 1 homolog-like; golgi SNARE 15 kDa protein; golgi SNARE with a size of 15 kDa; golgi integral membrane protein 3; vesicle transport protein GOS15; Bet1 golgi vesicular membrane trafficking protein like

Gene ID

[51272](#)

UniProt ID

[Q9NYM9](#)