

Product Information

MemDX™ Membrane Protein Human BRI3 (Brain protein I3) Full Length

Cat. No.: **MPC4584K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human BRI3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

BRI3

Protein Length

Full length

Protein Class

Receptor

TMD

2

Sequence

MDHKPLLQERPPAYNLEAGQGDIYACGPHGYGAIPAAAPPPYPYLVGTGIP
THHPRVYNIHSRTVTRYPAISIVVGGCPVCRVGVLEDCFTFLGIFLAI
LFPFGFICCFALRKRRCPNCGATFA

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

BRI3

Full Name

Brain protein I3

Alternative Names

BRI3; I3; pRGR2; Brain protein I3

Gene ID

[25798](#)

UniProt ID

[Q95415](#)