

Product Information

MemDX™ Membrane Protein Human COX16 (Cytochrome c oxidase assembly factor COX16)

Full Length

Cat. No.: **MPC4588K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human COX16 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

COX16

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MFAPAVMRAFRKNKTLGYGVPMLLLIVGGSFGLREFSQIRYDAVKSKMDP
ELEKKLKENKISLESEYEKIKDSKFDDWKNIRGPRPWEDPDLLQGRNPES
LTKTTT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

COX16

Full Name

Cytochrome c oxidase assembly factor COX16

Alternative Names

COX16; hCOX16; HSPC203; MC4DN22; C14orf112; cytochrome c oxidase assembly protein COX16 homolog, mitochondrial; COX16, cytochrome c oxidase assembly homolog; Cytochrome c oxidase assembly factor COX16

Gene ID

[51241](#)

UniProt ID

[Q9P0S2](#)