

Product Information

MemDX™ Membrane Protein Human RNF149 (Ring finger protein 149) Full Length

Cat. No.: **MPC4590K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human RNF149 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RNF149

Protein Length

Full length

Protein Class

Transferase

TMD

1

Sequence

MAWRRREASVGARGVLALALLALALCVPGARGRALEWFSAVVNIEYVDPQ
TNLTVWSVSESGRFGDSSPKEGAHGLVGVPWAPGGDLEGCAPDTRFFVPE
PGGRGAAPWVALVARGGCTFKDKVLVAARRNASAVVLYNEERYGNITLPM
SHAGTGNIVIMISYPKGREILELVQKGIPVTMTIGVGTRHVVQEFISGQS
VVFVAIAFITMMIISLAWLIFYIQRFLYTGSQIGSQSHRKETKKVIGQL
LLHTVKHGEKGIDVDAENCAVCIENTFKVKDIIRLPCKHIFHRICIDPWL
LDHRTCPMCKLDVIKALGYWGEPGDVQEMPAPESPPGRDPAANLSLALPD
DDGSDDSSPPSASPAESEPQCPSFKGDAGENTALLEAGRSDSRHGGPIS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

RNF149

Full Name

Ring finger protein 149

Alternative Names

RNF149; DNAPTP2; E3 ubiquitin-protein ligase RNF149; DNA polymerase-transactivated protein 2; RING-type E3 ubiquitin transferase RNF149; Ring finger protein 149

Gene ID

[284996](#)

UniProt ID

[Q8NC42](#)