

Product Information

MemDX™ Membrane Protein Human IGDCC3 (Immunoglobulin superfamily DCC subclass member 3) Full Length

Cat. No.: **MPC4595K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human IGDCC3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

IGDCC3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MAVQRAASPRRPPAPLWPRLLLPLLLLLLPAPSEGLGHS AELAFAVEPSD
DVAVPGQPIVLDCRVEGTPPV RITWRKNGVELPESTHSTLLANGSLMIRH
FRLEPGGSPSDEGDYECVAQNRFG LVVSRKARIQAATMSDFHVHPQATVG
EEGGVARFQCQIHGLPKPLITWEKNRVPIDTDNERYTLLPKGV LQITGLR
AEDGGIFHCVASNIASIRISHGARLTVSGSGSGAYKEPAILVGPENLT LT
VHQTAVLECVATGNPRPIVSWSR LDGRPIGVEGIQVLGTGNLIISDVTVQ
HSGVYVCAANRPGTRVRRTAQGR LVVQAPAEFVQHPQSISRPA GTTAMFT
CQAQGEPPPHVTWLKNGQVLGP GGHVRLKNNNSTLTISGIGPEDEAIYQC
VAENSAGSSQASARLTVLWAEGLPGPPRNRAVS SVSSTEVRVSWSEPLAN
TKEIIGYVLHIRKAADPPELEYQEAVSKSTFQHLVSDLEPSTAYSFYIKA
YTPRGASSASVPTLASTLGEAPAPPPLSVRVLGSSSLQLLWEPWPRLAQH
EGGFKLFYRPASKTSFTGPILLPGTVSSYNLSQLDPTAVYEVKLLAYNQH
GDGNATVRFVSLRGASERTALSPPCDCRKEEAANQTSTTGIVIGIHIGVT
CIIFCVLFLFFGQRGRVLLCKDVENQLSPPQGPRSQ RDPGILALNGARRG
QRGQLGRDEKRVDMKELEQLFP PASAAGQPDPRPTQDPAAPAPCEETQLS
VLPLQGCGLMGKTT EAKTTEATAPCAGLAAAPPPDGGPGLLSEGQASR
PAAARVTQPAHSEQ

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

IGDCC3

Full Name

Immunoglobulin superfamily DCC subclass member 3

Alternative Names

IGDCC3; PUNC; HsT18880; putative neuronal cell adhesion molecule; Immunoglobulin superfamily DCC subclass member 3

Gene ID

[9543](#)

UniProt ID

[Q8IVU1](#)