

Product Information

MemDX™ Membrane Protein Human RTP4 (Receptor transporter protein 4) Full Length

Cat. No.: **MPC4596K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human RTP4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RTP4

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MVVDFTWEQTFQELIQEAKPRATWTLKLDGNLQLDCLAQGWKQYQQRAF
GWFRCSQCQRSWASQVQILCHTYWEHWTSQGGQVRMRLFGQRCQKCSWSQ
YEMPEFSSDSTMRILSNLVQHILKKYYGNGTRKSPMPVILEVSLEGSHD
TANCEACTLGICGQGLKSCMTKPSKSLLPHLKTGNSSPGIGAVYLANQAK
NQSAAKEAKGSGYEKLGPSRDPDPLNICVFILLLVFIVVKCFTSE

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

RTP4

Full Name

Receptor transporter protein 4

Alternative Names

RTP4; IFRG28; Z3CXXC4; 28 kDa interferon-responsive protein; 28kD interferon responsive protein; 3CxxC-type zinc finger protein 4; receptor (chemosensory) transporter protein 4; zinc finger, 3CxxC-type 4; Receptor transporter protein 4

Gene ID

[64108](#)

UniProt ID

[Q96DX8](#)