

Product Information

MemDX™ Membrane Protein Human ILDR2 (Immunoglobulin like domain containing receptor 2) Full Length

Cat. No.: **MPC4598K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ILDR2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ILDR2

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MDRVLLRWISLFWLTAMVEGLQVTPDKKKVAMLFQPTVLRCHFSTSSHQ
PAVVQWKFKSYCQDRMGESLGMSSTRAQSLSKRNLEWDPYLDCLDSRRTV
RVVASKQGSTVTLGDFYRGREITIVHDADLQIGKLMWGDSGLYYCIITP
DDLEGNEDSVELLVLGRTGLLADLLPSFAVEIMPEWVFVGLVLLGVFLF
FVLVGICWCQCCPHSCCCYVRCPCPDSCCCPQALYEAGKAAKAGYPPSV
SGVPGPYSIPSVPLGGAPSSGMLMDKPHPPPLAPSDSTGGSHSVRKGYRI
QADKERDSMKVLYYVEKELAQFDPARRMRGRYNNITISELSSLHEEDSNFR
QSFHQMRSKQFPVSGDLESNPDYWSGVMGGSSGASRGPSAMEYNKEDRES
FRHSQPRSKSEMLSRKNFATGVPVSMDELAADFADSYGQRPRRADGNSHE
ARGGSRFERSESRAHSGFYQDDSL E EYQGQRSRSREPLTDADRGWAFSPA
RRRPAEDAHLPRLVSRTPGTAPKYDHSYLG SARERQARPEGASRGGSLET
PSKRSAQLGPRASYYAWSPPGTYKAGSSQDDQEDASDDALPPYSELELT
RGPSYRGRDLPYHSNSEKKRKKEPAKKTNDFPTRMSLVV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ILDR2

Full Name

Immunoglobulin like domain containing receptor 2

Alternative Names

ILDR2; C1orf32; dJ782G3.1; immunoglobulin-like domain-containing receptor 2; angulin-3; Immunoglobulin like domain containing receptor 2

Gene ID

[387597](#)

UniProt ID

[Q71H61](#)