

Product Information

MemDX™ Membrane Protein Human ZDHHC1 (Zinc finger DHHC-type containing 1) Full

Length

Cat. No.: **MPC4606K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC1

Protein Length

Full length

Protein Class

Transferase

TMD

4

Sequence

MYKMNICNKPSNKTAPEKSVWTAPAQPSGPSELQGQSRNRNGWSWPPHP
LQIVAWLLYLFFAVIGFGILVPLLPHHWVPAGYACMGAIFAGHLVVHLTA
VSIDPADANVRDKSYAGPLPIFNRSQHAHVIEDLHCNLCNVDVSARSKHC
SACNKCVCGFDHHCWLNNCVGERNYRFLHSVASALLGVLLLVLVATYV
FVEFFVNPMLRLRTNRHFEVLKNHTDVWFVFLPAAPVETQAPAILALAALL
ILLGLLSTALLGHLLCFHIYLMWHKLTTYEYIVQHRPPQEAKGVHRELES
CPPKMRPIQEMEFYMRTFRHMRPEPPGQAGPAAVNAKHSRPASPDPTPGR
RDCAGPPVQVEWDRKKPLPWRSPLLLLAMWGPQAPPCLCRKRGRGACIKC
ERLRPRIRRRGLGPPAAAPARRRIPRTPALCTPLALPAPTTRRRQSPWTR
FQWRRRAWAAPLWPPRGAGADSPRWRGRRVRPPFS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC1

Full Name

Zinc finger DHHC-type containing 1

Alternative Names

ZDHHC1; DHHC-1; ZNF377; C16orf1; HSU90653; palmitoyltransferase ZDHHC1; DHHC domain-containing cysteine-rich protein 1; probable palmitoyltransferase ZDHHC1; zinc finger DHHC domain-containing protein 1; zinc finger protein 377; zinc finger, DHHC domain containing 1; Zinc finger DHHC-type containing 1

Gene ID

[29800](#)

UniProt ID

[Q8WTX9](#)