

Product Information

MemDX™ Membrane Protein Human ABHD13 (Abhydrolase domain containing 13) Full

Length

Cat. No.: **MPC4611K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ABHD13 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ABHD13

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEKSWMLWNFVERWLIASWSWALCRISLLPLIVTFHLYGGIILLLLIF
ISIAGILYKFQDVLLYFPEQPSSSRLYVPMPTGIPHENIFIRTKDGIRLN
LILIRYTGDNPSYPTIYFHGNAGNIGHRLPNALLMLVNLKVNLLLV
RGYKSEGEASEEGLYLDSEAVLDYVMTRPDLDTKIFLFGRLGGAVAI
HLASENSHRISAIMVENTFLSIPHMASTLFSFFPMRYLPLWCYKNKFLSY
RKISQCRMPSLFISGLSDQLIPPVMMKQLYELSPSRTKRLAIFPDGTHND
TWQCQGYFTALEQFIKEVVKSHSPEEMAKTSSNVTII

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ABHD13

Full Name

Abhydrolase domain containing 13

Alternative Names

ABHD13; BEM46L1; C13orf6; bA153I24.2; protein ABHD13; abhydrolase domain-containing protein 13; alpha/beta hydrolase domain-containing protein 13; Abhydrolase domain containing 13

Gene ID

[84945](#)

UniProt ID

[Q7L211](#)