

Product Information

MemDX™ Membrane Protein Human AVL9 (AVL9 cell migration associated) Full Length

Cat. No.: **MPC4615K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human AVL9 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

AVL9

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEKARRGGDGVPRGPVLHIVVVGFFHHKGCQVEFSYPPLIPGDGHDSHTL
PEEWKYLPLALPDGAHNYQEDTVFFHLPPRNGNGATVFGISCYRQIEAK
ALKVRQADITRETVQKSVCLSKLPLYGLLQAKLQLITHAYFEEKDFSQI
SILKELYEHMNSSLGASLEGSQVYLGLSPRDLVHFRHKVLILFKLILL
EKKVLFYISPVNKLVGALMTVLSLFGMIEHGLSDCSQYRPRKSMSEDGG
LQESNPCADDFVSASTADVSHTNLGTIRKVMAGNHGEDAAMKTEEPLFQV
EDSSKGQEPNDTNQYLKPPSRPSPDSSSEDWETLDPVLEDPNLKEREQL
GSDQTNLFPKDSVPSESLPITVQPQANTGQVVLIPGLISGLEEDQYGMPL
AIFTKGYLCPLPYMALQQHHLLSDVTVRGFVAGATNILFRQQKHLSDAIVE
VEEALIQIHDPRLKLLNPTTADLRFADYLVRHVTENRDDVFLDGTGWEG
GDEWIRAQFAVYIHALLAATLQLDNEKILSDYGTTFVTAWKNTHNYRVWN
SNKHPALAEINPNHPFQGQYSVSDMKLRFHSVQNSERGGKIGNVMVTTTS
RNVVQTGKAVGQSVGGAFSSAKTAMSSWLSTFTTSTSQSLTEPPDEKP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

AVL9

Full Name

AVL9 cell migration associated

Alternative Names

AVL9; KIAA0241; late secretory pathway protein AVL9 homolog; AVL9 homolog (S. cerevisiae); AVL9 cell migration associated

Gene ID

[23080](#)

UniProt ID

[Q8NBF6](#)