

Product Information

MemDX™ Membrane Protein Human MFSD11 (Major facilitator superfamily domain containing 11) Full Length

Cat. No.: **MPC4622K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MFSD11 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MFSD11

Protein Length

Full length

Protein Class

Transporter

TMD

12

Sequence

MSPESKKLFNIIILGVAFMFMFTAFQTCGNVAQTVIRSLNRTDFHGSGYT
SMIIYGVFSASNLITPSVVAIVGPQLSMFASGLFYSMYIAVFIQPPFWS
FYTASVFIGIAAAVLWTAQGNCLTINSDEHSIGRNSGIFWALLQSSLFFG
NLYIYFAWQGKTQISESDRRTVFIALTVISLVGTVLFFLIRKPDSENVLG
EDESSDDQDMEVNESAQNNTKAVDAFKKSFKLCVTKEMLLSITTAYTG
LELTFFSGVYGTICGATNKFGAEEKSLIGLSGIFIGIGEILGGSFLGLLS
KNNRFGRNPVVLLGILVHFIAFYLI FLNMPGD APIAPVKGT DSSAYIKSS
KEVAILCSFLLGLGDSFNTQLLSILGFLYSEDSAPAF AIFKFVQSICAA
VAFFYSNYLLLHWQLLMVIFGFFGTISFFTVEWEEAAAFVARGSDYRSI

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MFSD11

Full Name

Major facilitator superfamily domain containing 11

Alternative Names

MFSD11; ET; UNC93-like protein MFSD11; major facilitator superfamily domain-containing protein 11; protein ET; Major facilitator superfamily domain containing 11

Gene ID

[79157](#)

UniProt ID

[O43934](#)