

Product Information

MemDX™ Membrane Protein Human RNF180 (Ring finger protein 180) Full Length

Cat. No.: **MPC4628K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human RNF180 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RNF180

Protein Length

Full length

Protein Class

Transferase

TMD

1

Sequence

MKRSKELITKNHSQEETSILRCWKCRKCIASSGCFMEYLENQVIKDKDDS
VDAQNICHVWHMNVEALPEWISCLIQKAQWTVGKLNCPFCGARLGGFNFV
STPKCSCGQLAAVHLSKSRDYQPTQAGRLMRPSVKYLSHPRVQSGCDKE
ALLTGGGSENRNHRLNLMARNNDPGRLTEALCLEVRPTYFEMKNEKLLS
KASEPKYQLFVPQLVTGRCATRAFHRSKSHSLDLNISEKLTLLPTLYEIH
KTTAYSRLNETQPIDLSGLPLQSSKNSYSFQNPSSFDPMSMLLQRFVAPH
ETQTQRGGEFQCGLEAASVYSDHTNTNLTFLMDLPSAGRSMPEASDQEE
HLSPLDFLHSANFSLGSINQRLNKRERSKLKLNLRKQRRRERWLQKQGY
SGVGLLDHMTLNNEMSTDEDNEYAEEKDSYICAVCLDVYFNPYMCYPCHH
IFCEPCLRTLAKDNPSSTPCPLCRTIISRVFFQTELNNATKTFFTKEYLK
IKQSFQKSNSAKWPLPSCRKAFLFGGFRRAAPVTRRQFPHGAHRMDYL
HFEDDSRGWWFDMDMVIIYISVNWWVIGFIVFCFLCYFFFPF

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

RNF180

Full Name

Ring finger protein 180

Alternative Names

RNF180; RINES; E3 ubiquitin-protein ligase RNF180; RING-type E3 ubiquitin transferase RNF180; Ring finger protein 180

Gene ID

[285671](#)

UniProt ID

[Q86T96](#)