

Product Information

MemDX™ Membrane Protein Human SHISA6 (Shisa family member 6) Full Length

Cat. No.: **MPC4629K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SHISA6 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SHISA6

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MALRRLLLLLLLLSLESLDLLPSVHGARGRAANRTLSAGGAAVGGRRAGGA
LARGGRELNGTARAPGIPEAGSRGQPAAAVAAAASAAVTYETCWGYDDV
SGQYDKEFECNNSSESGYLYCCGTCYYRFCCKKRHEKLDQRQCTNYQSPVW
VQTPSTKVVSPGPENKYDPEKDKNFTVYITCGVIAFVIVAGVFAKVSVD
KAHRPPREMNHRALADILRQQGPIPIAHCERETISAIDTSPKENTPVRS
SSKNHYTPVRTAKQTPEKPRMNNILTSATEPYDLSFSRSFQNLHLPPSY
ESAVKTNPSKYSSLKRLTDKEADEYYMRRRHLPDLAARGTLPLNVIQMSQ
QKPLPRERPRRPIRAMSQDRVLSADRGLPDEFSPYDRILSDEQLLSTER
LHSQDPLLSPERTAFPEQSLSRAISHTDVFVSTPVLDRYRMSKMHSHPSA
SNNSYATLGQSQTAAKRHAFASRRHNTVEQLHYIPGHHTCYTASKTEVTV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SHISA6

Full Name

Shisa family member 6

Alternative Names

SHISA6; protein shisa-6; protein shisa-6 homolog; shisa homolog 6; Shisa family member 6

Gene ID

[388336](#)

UniProt ID

[Q6ZSJ9](#)