

Product Information

MemDX™ Membrane Protein Human MGARP (Mitochondria localized glutamic acid rich protein) Full Length

Cat. No.: **MPC4640K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MGARP membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MGARP

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MYLRRAVSKTLALPLRAPNPAPLGKDASLRRMSSNRFPGSSGSNMIYYL
VVGVTVSAGGYAYKTVTSDQAKHTEHKTNLKEKTKAEIHPFQGEKENVA
ETEKASSEAPEELIVEAEVVDAAEESPSATVVVIKEASACPGHVEAAPETT
AVSAETGPEVTDAAARETTEVNPETTPEVTNAALDEAVTIDNDKDTTKNE
TSDEYAELEEENSPAEESESSAGDDLQEEASVGSEAASAQG

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MGARP

Full Name

Mitochondria localized glutamic acid rich protein

Alternative Names

MGARP; OSAP; HUMMR; CESP-1; C4orf49; protein MGARP; corneal endothelium-specific protein 1; hypoxia up-regulated mitochondrial movement regulator protein; ovary-specific acidic protein; Mitochondria localized glutamic acid rich protein

Gene ID

[84709](#)

UniProt ID

[Q8TDB4](#)