

Product Information

MemDX™ Membrane Protein Human MPV17L (MPV17 mitochondrial inner membrane protein like) Full Length

Cat. No.: **MPC4643K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MPV17L membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MPV17L

Protein Length

Full length

Protein Class

Receptor

TMD

4

Sequence

MAGWWPALSRAARRHPWPTNVLLYGSLVSAGDALQQRLQGREANWRQTRR
VATLVVTFHANFNYYVWLRLLEALPGRAPHALLAKLLCDQVVGAPIAVSA
FYVGMSILQGKDDIFDLKQKFWNTYLSGLMYWPFVQLTNFSLVPVQWRT
AYAGVCGFLWATFICFSQQSGDGTFSKSAFTILYTKGTSATEGYPKK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MPV17L

Full Name

MPV17 mitochondrial inner membrane protein like

Alternative Names

MPV17L; M-LPH; MLPH1; MLPH2; MPV17L1; mpv17-like protein; M-LP homolog; MPV17 mitochondrial membrane protein-like; Mpv17-like protein type 1; Mpv17-like protein type 2; MPV17 mitochondrial inner membrane protein like

Gene ID

[255027](#)

UniProt ID

[Q2QL34](#)