

Product Information

MemDX™ Membrane Protein Human MPV17L2 (MPV17 mitochondrial inner membrane protein like 2) Full Length

Cat. No.: **MPC4644K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MPV17L2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MPV17L2

Protein Length

Full length

Protein Class

Receptor

TMD

3

Sequence

MARGGWRRLLRRLSAGQLLFQGRALLVTNTLGCGALMAAGDGVQRQSWEIR
ARPGQVFDPRRSASMFVAGCSMGPFLLHYWYLSLDRLFPASGLRGFPNVLK
KVLVDQLVASPLLGVWYFLGLGCLEGGTVGESQCQLREKFWEFYKADWCV
WPAAQFVNFLFVPPQFRVTYINGLTGWDYLSYLYRSPVPLTPPGCVA
LDTRAD

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MPV17L2

Full Name

MPV17 mitochondrial inner membrane protein like 2

Alternative Names

MPV17L2; FKSG24; mpv17-like protein 2; MPV17 mitochondrial membrane protein-like 2; MPV17 mitochondrial inner membrane protein like 2

Gene ID

[84769](#)

UniProt ID

[Q567V2](#)