

Product Information

MemDX™ Membrane Protein Human SEL1L2 (SEL1L2 adaptor subunit of ERAD E3 ligase)

Full Length

Cat. No.: **MPC4646K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SEL1L2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SEL1L2

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MKPLSLLIEILILGVTIKTIKAEHNKRQKERNVTTQVSVNEIKQYLSH
ILEQRTSSNVINKRENLEKKKNQRKIRIKGIQNKDILKRKNHLQKQAE
KNFTDEGDQLFKMGIVLQQSKSQKQKEEAYLLFAKAADMGNLKAMEKMA
DALLFGNFGVQNITAAIQLYESLAKEGSCKAQNALGFLSSYGIGMEYDQA
KALIYYTFGSAGGNMMSQMLGYRYLSGINVLQNCEVALSYKKVADYIA
DTFEKSEGVPEKVRRLTERPENLSSNSEILDWDIYQYYKFLAERGDVQIQ
VSLGQLHLIGRKGLDQDYKALHYFLKAAKAGSANAMAFIGKMYLEGNA
VPQNNATAFKYFSMAASKGNAIGLHGLGLLYFHGKGVPLNYAEALKYFQK
AAEKGWPD AQFQLGFMYYSGSGIWKDYKLAFKYFYLASQSGQPLAIYYLA
KMYATGTGVVRSCTAVELYKGVCELGHWA EKFLTAYFAYKDG DIDSSLV
QYALLAEMGYEVAQSN S AFILESKKANILEKEKMYPMALLLWNRAAIQGN
AFARVKIGDYHYGYGTCKDYQTAATHYSIAANKYHNAQAMFNLAYMYEH
GLGITKDIHLARRLYDMAAQTSPDAHIPVLFVVMKLETT HLLRDILFFNF
TTRWNWLKLDNTIGPHWDLFVIGLIVPGLILLLRNHHG

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SEL1L2

Full Name

SEL1L2 adaptor subunit of ERAD E3 ligase

Alternative Names

SEL1L2; sel-1L2; C20orf50; protein sel-1 homolog 2; SEL1L2, ERAD E3 ligase adaptor subunit; sel-1 suppressor of lin-12-like 2; suppressor of lin-12-like protein 2; SEL1L2 adaptor subunit of ERAD E3 ligase

Gene ID

[80343](#)

UniProt ID

[Q5TEA6](#)