

Product Information

MemDX™ Membrane Protein Human KIAA2013 (KIAA2013) Full Length

Cat. No.: **MPC4661K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human KIAA2013 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

KIAA2013

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MWLQQRLKGLPGLSSSWARRLLCLLGLLLLLLWFGSGARRAAGGLHLL
PWSRGEPGAAEPSACLEAATRAWRLRERGEVVPLGPGVPALVANGFLAL
DVAANRLWVTPGEREPAVAPDFVPFVQLRPLSALAEAGEAVLLLREGLLR
RVRCLQLGSPGPGPVAAGPGPASVSGLAAGSGRDCVLLQEDFLAHRGRPH
VYLQRIQLNPNPTERVAALQTVGPTAGPAPKAFTSTLEKVGDPHQFLYSGR
SPPTPTGLVHLVVAAKKLVNRLQVAPKTQLDETVLWVHVSGPINPQVL
KSKAAKELKALQDLARKEMLELLDMPAAELLQDHQLLWAQLFSPGVEMKK
ITDHTPSGLTVNLTYYMLSCSPAPLLSPSLSHRERDQMESTLNYEDHC
FSGHATMHAENLWPGRLLSSVQQILQLSDLWRLTLQKRGCKGLVKVGAPGI
LQGMVLSFGGLQFTENHLQFQADPDVLHNSYALHGIRYKNDHINLAVLAD
AEGKPYLHVSVESRGQPVKIYACKAGCLDEPVELTSAPTGHTFSVMVTQP
ITPLLYISTDLTHLQDLRHTLHLKAILAHDEHMAQQDPGLPFLWFSVAS
LITLFLFLFKLIYNEYCGPGAKPLFRSKEDPSV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

KIAA2013

Full Name

KIAA2013

Alternative Names

KIAA2013; uncharacterized protein KIAA2013; RP5-1077B9.1; KIAA2013

Gene ID

[90231](#)

UniProt ID

[Q8IYS2](#)