

Product Information

MemDX™ Membrane Protein Human TARM1 (T cell-interacting, activating receptor on myeloid cells 1) Full Length

Cat. No.: **MPC4664K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human TARM1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TARM1

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MIPKLLSLLCFRLCVGQGDTRGDGLPKPSLSAWPSSVVPANSNVTLRCW
TPARGVSFVLRKGGIILESPKPLDSTEGAAEFHLNNLKVRNAGEYTCYY
RKASPHILSQRSDVLLLLVTGHLSPFLRTYQRGTVTAGGRVTLQCQKRD
QLFVPIMFALLKAGTPSPIQLQSPAGKEIDFSLVDVTAGDAGNYSCMYQ
TKSPFWASEPSDQLEILVTVPPTSSNYSLGNFVRLGLAAVIVVIMGAF
LVEAWYSRNVSPGESEAFKPE

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TARM1

Full Name

T cell-interacting, activating receptor on myeloid cells 1

Alternative Names

TARM1; OLT-2; T-cell-interacting, activating receptor on myeloid cells protein 1; OSCAR-like transcript-2 protein; T cell-interacting, activating receptor on myeloid cells 1

Gene ID

[441864](#)

UniProt ID

[B6A8C7](#)