

Product Information

MemDX™ Membrane Protein Human C4orf3 (Chromosome 4 open reading frame 3) Full

Length

Cat. No.: **MPC4676K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human C4orf3 membrane protein expressed in *E.coli*. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

C4orf3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEVDAPGVDGRDGLRERRGFSEGGQRQNFVDPQSGANGLPKHSYWLDLWL
FILFDVVVFLFVYFLP

Product Description

Expression Systems

E.coli

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

C4orf3

Full Name

Chromosome 4 open reading frame 3

Alternative Names

C4orf3; ALN; HCVFTP1; uncharacterized protein C4orf3; HCV F-transactivated protein 1; another-regulin; hepatitis C virus F protein-transactivated protein 1; Chromosome 4 open reading frame 3

Gene ID

[401152](#)

UniProt ID

[Q8WVX3](#)