

Product Information

MemDX™ Membrane Protein Human NSG2 (Neuronal vesicle trafficking associated 2) Full Length

Cat. No.: **MPC4684K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human NSG2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

NSG2

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MVKLNSNPSEKGTKPPSVEDGFQTVPLITPLEVNHQLQPAPEKVIVKTRT
EYQPEQKNKGKFRVPKIAEFTVTILVSLALAFLACIVFLVVYKAFTYDHS
CPEGFVYKHKRCIPASLDAYYSSQDPNSRSFYTVISHYSVAKQSTARAI
GPWLSAAAVIHEPKPPKTQGH

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

NSG2

Full Name

Neuronal vesicle trafficking associated 2

Alternative Names

NSG2; CALY3; HMP19; neuronal vesicle trafficking-associated protein 2; calcyon family member 3; hypothalamus golgi apparatus expressed 19 kDa protein; neuron-specific protein family member 2; p19 protein; protein p19; Neuronal vesicle trafficking associated 2

Gene ID

[51617](#)

UniProt ID

[Q9Y328](#)