

Product Information

MemDX™ Membrane Protein Human NPIP3 (Nuclear pore complex interacting protein family member B3) Full Length

Cat. No.: **MPC4688K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human NPIP3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

NPIP3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MVKLSIVLTPQFLSHDQGQLTKELQQHVKSVCPCYLRKVINTLADHHH
RGTDGFGSPWLHVIAFPTS YKVVITLWIVYLWVSLKTIFWSRNGHDGS
TDVQQRRAWRSNRRRQEGLRSICMHTKKRVSSFRGNKIGLKDVITLRRHVE
TKVRAKIRKRKVTTKINHHDKINGKRKTARKQKMFQRAQELRRRAEDYHK
CKIPPSARKALCNWVRMAAAEHRHSSGLPYWPLYTAETLKNRMGHQPPPP
TQQHSITDNSLSLKTPECLLTPLPPSADDNLKTPPECVLTPLPPSADDN
LKTPPECVLTPLPPSADDNLKTPPECCLTPLPPSADDNLKTPPECCLTPL
PPSALPSAPPSADDNLKTRAECCLHPLPPSADDNLKTPSERQLTPLPPSA
PPSADDNIKTPAERLRGPLPPSADDNLKTPSERQLTPLPPSAPPSADDNI
KTPAERLRGPLPPSADDNLKTPSERQLTPLPPSAPPSADDNIKTPAERLR
GPLPPSADDNLKTPSERQLTALPPSADDNIKTPAERLRGPLPPSADDNLK
TPSERQLTPLPPSAPPSADDNIKTPAFHPQRMII SRHLP SVSSLPFHPQL
HPQQMIISRHLPSVCGGRFHPQQMIISRHLPSVSSLPFHPQLHPQQMIIS
RHLPSVCGGRFHPQRMII SRHLP SVSSLPFHPQLHPQQMIISRHLPSVCG
GRFHPQQMIISRHLPSVSSLPFHPQLHPQQMIISRHLPSVCGGRFHPQRM
IISRHLPSVSSLPFHPQLHPQQMIISRHLPSVCGERLRGPLPPSADDNLK
TPSERQLTPLPPSAPPSADDNIKTPAERLRGPLPPSADDNLKTPSERQLT
PLPPSAPPSADDNIKTPAERLRGPLPPSADDNLKTPSERQLTPLPPSAPP
SADDNIKTPAERLRGPLPPSADDNLKTPPLATQEA EAEKPRKPKRQRAAE
MEPPPEPKRRRVGDVEPSRKPKRRRAADVEPSSPEPKRRRVGDVEPSRK
KRRRAADVEPSSPEPKRRRVGDVEPSRKPKRRRAADVEPSLPEPKRRRLS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

NPIP3

Full Name

Nuclear pore complex interacting protein family member B3

Alternative Names

NPIP3; NPIP; NPIPB; NPIPB5; NPIPL3; KIAA0220L; KIAA0220-like protein; nuclear pore complex-interacting protein B type; nuclear pore complex-interacting protein-like 3; protein pps22-1; Nuclear pore complex interacting protein family member B3

Gene ID

[23117](#)

UniProt ID

[Q92617](#)