

# Product Information

## MemDX™ Membrane Protein Human FAM234A (Family with sequence similarity 234 member

### A) Full Length

Cat. No.: **MPC4689K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human FAM234A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

FAM234A

#### Protein Length

Full length

#### Protein Class

Receptor

#### TMD

1

#### Sequence

MLDHKDLEAEIHPLKNEERKSQENLGNPSKNEDNVKSAPPQSRLSRCRAA  
AFFLSLFLCLFVVFVVSFVIPCPDRPASQRMWRIDYSAAVIYDFLAVDDI  
NGDRIQDVLFLYKNTNSSNNFSRSCVDEGFSSPCTFAAAVSGANGSTLWE  
RPVAQDVALVECAVPQPRGSEAPSACILVGRPSSFIAVNLFTGETLWNHS  
SSFSGNASILSPLLQVPDVGDDGAPDLLVLTQEREEVSGHLYSGSTGHQI  
GLRGS LGVDGESGFLHVT RTGAHYILFPCASSLCGCSVKGLYEKVTGSG  
GPFKSDPHWESMLNATTRMLSHSSGAVRYLMHVPGNAGADVLLVGSEAF  
VLLDGQELTPRWTPKAAHVLRKPIFGRYKPD TLAVAVENG TGTD RQILFL  
DLGTGAVLCSLALPSLPGGPLSASLPTADHRSAFFFWGLHELGSTSETET  
GEARHSLYMFHPTLPRVLELANVSTHIVAFDAVLFEPSRHAAYILLTGP  
ADSEAPGLVSVIKHKVRDLVPSSRVVRLGEGGPDSQAIRDRFSRLRYQS  
EA

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

**Protein Format**

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

FAM234A

**Full Name**

Family with sequence similarity 234 member A

**Alternative Names**

FAM234A; gs19; ITFG3; C16orf9; protein FAM234A; integrin alpha FG-GAP repeat containing 3; protein ITFG3; Family with sequence similarity 234 member A

**Gene ID**

[83986](#)

**UniProt ID**

[Q9H0X4](#)