

Product Information

MemDX™ Membrane Protein Human ZDHHC11B (Zinc finger DHHC-type containing 11B)

Full Length

Cat. No.: **MPC4706K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC11B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC11B

Protein Length

Full length

Protein Class

Transferase

TMD

5

Sequence

MDTRSGSQCSVTPEAIRNNEELVLPPRISRVNGWSLPLHYFRVVTWAVFV
GLSLATFRIFIPLLPHSWKYIAYVVTGGIFSFLVHVLHASCIDPADSNV
RLMKNYSQPMPLFDRSKHAHVQNQFCHLCKVTVNKTKHCISCNKCMSG
FDHHCKWINNCVGSRNWFFFSTVASATAGMLCLIAILLYVLVQYLVNPR
VLRTDPRYEDVKNMNTWLLFLPLFPVQVQTLIVVIIRMLVLLDLLGLVQ
LGQLLIFHIYKAKKMTTFEYLINTRKEESSKHQAVRKDPYVQMDKGFLQ
QGAGALGSSAQGVKAKSSLLIYKCPCHFCTSVNQDGDGSKAQEADDAPSTS
TLGLQQETTEPMKTDSESED

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC11B

Full Name

Zinc finger DHHC-type containing 11B

Alternative Names

ZDHHC11B; DHHC-11B; probable palmitoyltransferase ZDHHC11B; zinc finger DHHC domain-containing protein 11B; Zinc finger DHHC-type containing 11B

Gene ID

[653082](#)

UniProt ID

[P0C7U3](#)