

Product Information

MemDX™ Membrane Protein Human RNF183 (Ring finger protein 183) Full Length

Cat. No.: **MPC4711K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human RNF183 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RNF183

Protein Length

Full length

Protein Class

Transferase

TMD

1

Sequence

MAEQQGREGLEAECPCWNPFNNTFHTPKMLDCCHSFCVECLAHLSLVTPA
RRRLLCPLCRQPTVLASGQPVTDLPTDTAMLALLRLEPHHVILEGHQLCL
KDQPKSRYFLRQPQVYTLDLGPQPGGQTGPPPDATVSTPILIPSHHS
LRECFRNPQFRIFAYLMAVILSVTLLIFSIWTKQFLWGVG

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

RNF183

Full Name

Ring finger protein 183

Alternative Names

RNF183; E3 ubiquitin-protein ligase RNF183; probable E3 ubiquitin-protein ligase RNF183; Ring finger protein 183

Gene ID

[138065](#)

UniProt ID

[Q96D59](#)