

Product Information

MemDX™ Membrane Protein Human RNF182 (Ring finger protein 182) Full Length

Cat. No.: **MPC4712K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human RNF182 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RNF182

Protein Length

Full length

Protein Class

Transferase

TMD

2

Sequence

MASQPPEDTAESQASDELECKICYNRYNLKQRKPKVLECCHRVCAKCLYK
IIDFGDSPQGIVCPFCRFETCLPDDEVSSLPDDNNILVNLTCGGKGKKC
LPENPTELLTPKRLASLVSPSHTSSNCLVITIMEVQRESSPSLSSTPVV
EFYRPASFDSVTTVSHNWTWVNCTSLLFQTSIRVLVWLLGLLYFSSLPLG
IYLLVSKKVTLGVVVSVSLVPSSLVILMVYGFCQCVCHEFLDCMAPPS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

RNF182

Full Name

Ring finger protein 182

Alternative Names

RNF182; E3 ubiquitin-protein ligase RNF182; RING-type E3 ubiquitin transferase RNF182; Ring finger protein 182

Gene ID

[221687](#)

UniProt ID

[Q8N6D2](#)