

Product Information

MemDX™ Membrane Protein Human CTXN3 (Cortixin 3) Full Length

Cat. No.: **MPC4745K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CTXN3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CTXN3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MDGGQPIPSLVPLGNESADSSMSLEQKMTFVFVILLFIFLGILIVRCFR
ILLDPYRSMPTSTWADGLEGLEKGQFDHALA

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw

cycles.

Target

Target Protein

CTXN3

Full Name

Cortixin 3

Alternative Names

CTXN3; KABE; cortixin-3; kidney and brain-expressed protein; Cortixin 3

Gene ID

[613212](#)

UniProt ID

[Q4LDR2](#)