

Product Information

MemDX™ Membrane Protein Human PAPPA-AS1 (PAPPA antisense RNA 1) Full Length

Cat. No.: **MPC4756K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human PAPPA-AS1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PAPPA-AS1

Protein Length

Full length

Protein Class

Receptor

TMD

2

Sequence

MRYGFVRKKHRGLFLTVAALPIWNPISFVKWYKSHKLSQHCIRICGHL
CQKHLDMFLSVIGQRWPIDVFSSVFDHQVSAIGSDIWWFLKLFLVSFFF
FF

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

PAPPA-AS1

Full Name

PAPPA antisense RNA 1

Alternative Names

PAPPA-AS1; DIPAS; TAF4B; PAPPAS; TAF2C2; PAPPAAS; PAPPA-AS; TAFII105; NCRNA00156; DIPLA1-antisense expressed; PAPPA antisense RNA 1 (non-protein coding); PAPPA antisense RNA 1 (tail to tail); PAPPA antisense RNA 1

Gene ID

[493913](#)

UniProt ID

[Q5QFB9](#)