

Product Information

MemDX™ Membrane Protein Human IFI27L2 (Interferon alpha inducible protein 27 like 2)

Full Length

Cat. No.: **MPC4759K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human IFI27L2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

IFI27L2

Protein Length

Full length

Protein Class

Receptor

TMD

3

Sequence

MMKRAAAAVGGALAVGAVPVVLSAMGFTGAGIAASSIAAKMMSAAAIAN
GGGVSAGSLVATLQSVGAAGLSTSSNILLASVGSVLGACLGNSPSSSLPA
EPEAKEDEARENVPQGEPPKPPLKSEKHEE

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

IFI27L2

Full Name

Interferon alpha inducible protein 27 like 2

Alternative Names

IFI27L2; TLH29; FAM14A; ISG12B; interferon alpha-inducible protein 27-like protein 2; ISG12(b) protein; family with sequence similarity 14, member A; interferon-stimulated gene 12b protein; pIFI27-like protein; Interferon alpha inducible protein 27 like 2

Gene ID

[83982](#)

UniProt ID

[Q9H2X8](#)