

Product Information

MemDX™ Membrane Protein Human LINC00052 (Long intergenic non-protein coding RNA

52) Full Length

Cat. No.: **MPC4777K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human LINC00052 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

LINC00052

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MNCDALLHHS AIPEDFLHIFLLQKISVSLPLSLSQSVCLFYISISLCVSL
LLHISLCVSVYVSLSLSSFPCLTHTHTHSQSKDTSVLTFTFCFKQHT
HFTLNYTSHAHELSAPSVHPTCVFTFKAAPSPRPAT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

LINC00052

Full Name

Long intergenic non-protein coding RNA 52

Alternative Names

LINC00052; TMEM83; NCRNA00052; transmembrane protein 83; Long intergenic non-protein coding RNA 52

Gene ID

[145978](#)

UniProt ID

[Q96N35](#)