

Product Information

MemDX™ Membrane Protein Human PXMP4 (Peroxisomal membrane protein 4) Full Length

Cat. No.: **MPC4814K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human PXMP4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PXMP4

Protein Length

Full length

Protein Class

Receptor

TMD

2

Sequence

MAAPPQLQALLQAVNKLLRQRRYHAALAVIKGFRNGAVYGVKIRAPHALV
MTFLFRSGSLREKLQAILKATYIHSRNLACFVFAYKSLHALQSHVQGETH
QMHSFLAAFIFGGLLLFGENNNINSQINMYLTSRVLYALCRLGVEKGYIPA
LKWDPFPLHTAVIWGLVLWLFHEYHRPTLQPSLQSSMTYLYEDSNVWHDLS
DFLIFNKSHPSK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

PXMP4

Full Name

Peroxisomal membrane protein 4

Alternative Names

PXMP4; PMP24; 24 kDa peroxisomal intrinsic membrane protein; peroxisomal membrane protein 4, 24kDa; Peroxisomal membrane protein 4

Gene ID

[11264](#)

UniProt ID

[Q9Y6I8](#)