

# Product Information

## MemDX™ Membrane Protein Human CLEC2L (C-type lectin domain family 2 member L) Full Length

Cat. No.: **MPC4818K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CLEC2L membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

CLEC2L

#### Protein Length

Full length

#### Protein Class

Receptor

#### TMD

1

#### Sequence

MEPAREPPSRARPPPPPLAARPAPAPAAPRPRSPAEAEARGPEGLLRSGS  
GYEGSTSWKAALEDTTTLLLGAIIVLLFAILVVMISILASKGCIKCEAPC  
PEDWLLYGRKCYFFSEEPDWNTGRQYCHTHEAVLAVIQSQKELEFMFKF  
TRREPWIGLRRVGDEFHWVNGDPFDPDTFTIAGPGECVFVEPTRLVSTEC  
LMTRPWVCSKMAYT

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

#### Form

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

CLEC2L

**Full Name**

C-type lectin domain family 2 member L

**Alternative Names**

CLEC2L; FLJ32986; C-type lectin domain family 2 member L

**Gene ID**

[154790](#)

**UniProt ID**

[P0C7M8](#)