

Product Information

MemDX™ Membrane Protein Human TSPAN19 (Tetraspanin 19) Full Length

Cat. No.: **MPC4834K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human TSPAN19 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TSPAN19

Protein Length

Full length

Protein Class

Receptor

TMD

3

Sequence

MLRNNKTIHKYFLNLINGAFLVLGLLFMGFGAWLLDRNNFLTAFDENN
HFIVPISQILIGMGSSTVLFCLLGYGIIHNEIRWLLIVYAVLITWTFVAVQ
VVLSAFIITKKEEVQQLWHDKIDFVISEYGSKDKPEDITKWITLNALQKT
LQCCGQHNYTDWIKNKNKENSQVPCSTKSTLRKWFCDEPLNATYLEGC
ENKISAWYNVNVLTIGINFGLLTSEVFQVSLTVCFFKNIKNIIHAEM

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TSPAN19

Full Name

Tetraspanin 19

Alternative Names

TSPAN19; tetraspanin-19; putative tetraspanin-19; tspan-19; Tetraspanin 19

Gene ID

[144448](#)

UniProt ID

[P0C672](#)