

Product Information

MemDX™ Membrane Protein Human C17orf78 (Chromosome 17 open reading frame 78) Full Length

Cat. No.: **MPC4843K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human C17orf78 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

C17orf78

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MDTILVFSLIIASYDANKKDLRDSSCRLEQLPGIFPKDVRSIRELQMQET
HTETKRITTFIQNRTIATLQCLGSDSKVKVNLVYLERRPKVKHILKNLRRI
AAPRRNSSASSSCHLIPTSKFQTGSLKKGKAFLPGISQCKVLGASSETFP
TTAPSITPGNKEGEKTTSTDTDENLEKRQKWSIVVKILIAVTLTLLSGVAI
IVFVIFEVPCPYQCLGARKLCQCQWLWRWQKKGGQPPGTAESKPDSQPQK
VGQDAANSSNPKKAAEITVIHQTYF

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

C17orf78

Full Name

Chromosome 17 open reading frame 78

Alternative Names

C17orf78; uncharacterized protein C17orf78; Chromosome 17 open reading frame 78

Gene ID

[284099](#)

UniProt ID

[Q8N4C9](#)