

# Product Information

## MemDX™ Membrane Protein Human FAM118A (Family with sequence similarity 118 member

### A) Full Length

Cat. No.: **MPC4860K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human FAM118A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

FAM118A

#### Protein Length

Full length

#### Protein Class

Receptor

#### TMD

1

#### Sequence

MDSVEKTTNRSEQKSRKFLKSLIRKQPQELLLVIGTGVSAAVAPGIPALC  
SWRSCIEAVIEAAEQLEVLHPGDVAEFRRKVTKDRDLLVVAHDLIRKMSP  
RTGDAKPSFFQDCLMEVFDDLEQHIRSPVVLQSILSLMDRGAMVLTNNYD  
NLLEAFGRRQNKPMESLDLKDGTKVLEWARGHMKYGVLIHGLYTDPCGV  
VLDPSGYKDVTQDAEVMEVLQNLRYRTKSFLVGCGETLRDQIFQALFLYS  
VPNKVDLEHYMLVLKENEDHFFKHQADMLLHGKVVSYGDCFDHFPGYVQ  
DLATQICKQQSPDADRVDSTLLGNACQDCAKRKLEENGIEVSKKRTQSD  
TDDAGGS

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

FAM118A

**Full Name**

Family with sequence similarity 118 member A

**Alternative Names**

FAM118A; C22orf8; protein FAM118A; bK268H5.C22.4; Family with sequence similarity 118 member A

**Gene ID**

[55007](#)

**UniProt ID**

[Q9NWS6](#)