

Product Information

MemDX™ Membrane Protein Human FRMD5 (FERM domain containing 5) Full Length

Cat. No.: **MPC4898K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human FRMD5 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FRMD5

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MLSRLMSGSSRSLEREYSCTVRLDDSEYTCTIQRDAKGQYLFDLLCHHL
NLLEKDYFGIRFVDPDKQRHWLEFTKSVVKQLRSQPPFTMCFRVKFYPAD
PAALKEEITRYLVFLQIKRDLYHGRLCKTSDAALLAAYILQAEIGDYDS
GKHPEGYSSKFQFFPKHSEKLERKIAEIHKTELSGQTPATSELNFLRKAQ
TLETYGVDPHPCCKDVSGNAFLAFTPFGFVVLQGNKRVHFIKWNEVTKLK
FEGKTFYLYVSQKEEKIILTYFAPTPEACKHLWKCGIENQAFYKLEKSS
QVRTVSSSNLFFKGSFRFRYSGRVAKEVMESSAKIKREPPEIHRAGMVPSR
SCPSITHGPRLSSVPRTRRRRAVHISIMEGLESLRDSAHSTPVRSTSHGDT
FLPHVRSSRTDSNERVAVIADEAYSPADSVLPTPVAEHSLELMLLSRQIN
GATCSIEEEKESEASTPTATEVEALGGELRALCQGHSGPEEEQVNKFVLS
VLRLLLVTMGLLFVLLLLLILTESDLIAFFRDIRQTPEFEQFHYQYFC
PLRRWFACKIRSVVSLIDT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

FRMD5

Full Name

FERM domain containing 5

Alternative Names

FRMD5; FERM domain-containing protein 5; FERM domain containing 5

Gene ID

[84978](#)

UniProt ID

[Q7Z6J6](#)