

Product Information

NONIDET P40 Substitute

Igepal CA-630; [Octylphenoxy]Polyethoxyethanol; NONIDET P40 Substitute

Cat. No.: **MPD0169K**

This product is for research use only and is not intended for diagnostic use.

MPD0169K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Detergent

Classification

MB Reagents

Molecular Formula

C₃₂H₅₈O₁₀

CAS Number

2497-59-8

Formula Weight

603

CMC in H₂O

0.05-0.3 mM

Size

500 ml; 1 L

Storage

Aliquot and store at -20°C.