

## Product Information

### Pluronic F-68

**Polyoxyethylene-Polyoxypropylene Block Copolymer; Methyl-Oxirane, Polymer with Oxirane, (C<sub>3</sub>H<sub>6</sub>O.C<sub>2</sub>H<sub>4</sub>O)<sub>x</sub>; Poloxamer 188; Pluronic F-68**

Cat. No.: **MPD0175K**

This product is for research use only and is not intended for diagnostic use.

MPD0175K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Type

Detergent

### Classification

MB Reagents

### Molecular Formula

EO78PO30EO78

### CAS Number

9003-11-6

### Formula Weight

8400

### CMC in H<sub>2</sub>O

17.9 mM

### Size

100g; 500g

### Storage

Aliquot and store at -20°C.