

Product Information

Fos-Choline-14, Deuterated

n-Tetradecylphosphocholine-d42; Fos-Choline-14, Deuterated

Cat. No.: **MPD0193K**

This product is for research use only and is not intended for diagnostic use.

MPD0193K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Detergent

Classification

Deuterated

Molecular Formula

C₁₉D₄₂NO₄P

CAS Number

869638-98-2

Formula Weight

421.5

CMC in H₂O

5.1E-05

Purity

≥ 90% by HPLC analysis

Solubility in water

≥ 1% (in water at 20°C)

Size

100mg; 500mg; 1g

Storage

Aliquot and store at -20°C.