

Product Information

Fos-Choline-8, Fluorinated

(1H, 1H, 2H, 2H-Perfluorooctyl)phosphocholine; Fos-Choline-8, Fluorinated

Cat. No.: **MPD0198K**

This product is for research use only and is not intended for diagnostic use.

MPD0198K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Detergent

Classification

Fluorinated

Molecular Formula

C₁₃H₁₇F₁₃NO₄P

CAS Number

313997-22-7

Formula Weight

529.2

CMC in H₂O

2.9 mM

Purity

≥ 99% by HPLC analysis

Solubility in water

≥ 10% (in water at 20°C)

Conductance

< 200 μS

Size

1g; 5g; 25g

Storage

Aliquot and store at -20°C.