

Product Information

L-(+)-Selenomethionine

(S)-2-Amino-4-(Methylseleno)butyric Acid; L-(+)-Selenomethionine

Cat. No.: **MPD0223K**

This product is for research use only and is not intended for diagnostic use.

MPD0223K-1.png

This product is a detergent for membrane protein solubilization. It is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Type

Detergent

Classification

Selenium

Molecular Formula

C₅H₁₁NO₂Se

CAS Number

3211-76-5

Formula Weight

196.1

Purity

≥ 98% by HPLC analysis

Size

250mg; 1g; 5g; 25g

Storage

Aliquot and store at -20°C.